

DETAIL & FEATURES PILOT SERIES



LCD Touch Screen

LCD Touch Screen touch panel allows easy and convenient control. With data analysis by graphic trend for automatic monitoring.



25
1978.0 m

Temperature recorder

Current status is being monitored. Any event must be recorded and reported, program status is displayed.



Drying Chamber with condenser

Sophisticated design for fine quality drying results. *Stoppering(Option): Sealing is completed under vacuum condition (0.3mBar) automatically.

Lyoph-pride 03					
External Dim.(mm)	(W)700 x (D)750 x (H)865	Shelf Control range	-45 ~ + 60°C		
Internal Dim.(mm)	(W)360 x (D)470 x (H)228	Condenser	-70℃		
Shelf Dim.(mm)	(W)275 x (D)345 x (H)16	Defrost	Hotgas control		
Shelf Area	2+1,0.09m² (EA)	#15hinBioBranch 527271 · 초저온냉동고 · 제빙기 · 혈액냉동고			
Clearance	70mm				
Collect Cap.(Liters)	4L	_{팀장(차장)} 김 태 성 (보	╒) 일신 바이오브랜치		
Option	Stoppering		2118 서울특별시 중랑구 중랑천로 77 통전화 : (02) 491-1327, 1328, 1329 스 : (02) 491-4073		
Οριοπ	manifold chamber	© 1577-4053 E-	-mail:isbb4053@hanmail.net obile:010-4577-9906		

 Ising
 (주) 일신바이오베이스
 대표번호
 1577-4053

 (주) 일신바이오브랜지
 직통번호
 02-491-1327~1329